FOLLOW-UP

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EXPAND-GRAFT® PTFE VASCULAR GRAFT

Description of the product: The EXPAND-GRAFT® PTFE vascular graft consists of a white, cylindrical, flexible tube, with a printed guideline to direct its implantation. The integral water permeability is $\leq 1 \text{ ml/cm}^2/\text{min}$. EXPAND-GRAFT® is supplied sterile, individually packed in double blister and placed in a cardboard box. They are available in a wide range of sizes, in straight or tapered shapes. The following sizes are available (D = diameter, L = Length)

D x L	10 cm	20 cm	25 cm	30 cm	40 cm	50 cm	60 cm	70 cm	80 cm
3,5 mm	20310	20320	20325	20330	20340	20350	20360	20370	20380
4 mm	20410	20420	20425	20430	20440	20450	20460	20470	20480
5 mm	20510	20520	20525	20530	20540	20550	20560	20570	20580
6 mm	20610	20620	20625	20630	20640	20650	20660	20670	20680
7 mm	20710	20720	20725	20730	20740	20750	20760	20770	20780
8 mm	20810	20820	20825	20830	20840	20850	20860	20870	20880
9 mm	20910	20920	20925	20930	20940	20950	20960	20970	20980
10 mm	21010	21020	21025	21030	21040	21050	21060	21070	21080
4/6 mm	-	-	-	-	24640	24650	-	24670	-
4/7 mm	-	-	-	-	24740	24750	-	24770	-
6/8 mm	-	-	-	-	26840	26850	-	26870	-

Indications:

Repair or substitution of peripheral arteries.
Extraanatomical bypass
Arteriovenous vascular access

Content: 1 product unit, instructions for use leaflet, follow-up form and traceability labels.

Contraindications:

EXPAND-GRAFT® is not indicated for use in the coronary arteries
Vascular grafts should not be used in patients with anatomical, physiological or pathological conditions which can prevent the use of the material, or that may result in high risk during the procedure, such as: cardiac, respiratory or renal insufficiency, recent myocardial infarction, obesity, advanced age, extreme frailty, cardiopathies, neoplasia or low life expectancy.

Warnings and precautions:

- Product should be stored in a dry place, between 0 and 50°C.
- Keep away from sunlight.
- Do not peel, skin, separate or delaminate any portion of the graft wall.
- Do not expose the graft to organic solvent solutions such as alcohol, and do not force aqueous solutions through the graft wall

- The graft should always be handled with clean, sterile gloves.
- Insert suture at least 2 mm away from the end of the graft
- Excessive tension on the graft can lead to anastomotic rupture, however, moderate tension is desired.
- In vascular access procedures, do not insert the cannula less than 2.5 cm away from the anastomosis.
- Puncture sites should be rotated and spaced, since multiple punctures in the same area can damage the graft.
- Punctures should be inserted in a 45° angle, to avoid damaging the opposite wall.
- Do not twist the graft, use the guideline as reference.
- Always use sharp instruments to cut the graft.
- Preferably use monofilament non absorbable sutures with taper point.
- Always use atraumatic or covered clamps or clips.
- Do not use hemostatic clips in vascular access procedures.
- Do not use if package is damaged.
- Single use product. Do not reuse.
- Sterile. Do not resterilize. Sterilized by ETO.

Adverse reactions: Adverse reactions and potential complications associated with any vascular graft include, but are not limited to: occlusion, coagulopathy, infection, thromboembolism, aneurysm, neointimal hyperplasia, bleeding, seroma weeping. These complications, together or separately can lead to amputation or death.

Possible symptoms related to these adverse reactions include, but are not limited to: fever, shivering, pain, blushing, swelling, bleeding, haematoma in the operated limb, numbness, low temperature of the operated limb, chest pain or pressure, nausea or vomit, excessive sweating, dizziness and fainting. The occurrence of any of these symptoms should be communicated to the physician, and the patency of the graft should be assessed periodically. Replacement or clearing of the graft may be necessary.

Procedures for post-surgical clearance:

Embolectomy or angioplasty balloon catheter usage: The size of the balloon should be carefully chosen, since the graft does not distend. Excessively large catheters can damage the graft or rupture the balloon.

Thrombectomy: Thrombectomy procedures can be performed by carrying out a longitudinal incision with stay sutures. When a transverse incision is used, a horizontal mattress suture technique may be required to aid closure.

When inserting a catheter in the graft, its diameter should never be larger than the grafts lumen to avoid damage.

Traceability: Traceability labels are supplied attached to the instructions for use leaflet. The follow-up form at the end of this leaflet is the official channel of communication with the manufacturer regarding product performance. It should be filled entirely, describing the occurrence and physician's feedback, with one of the traceability labels stucked to it in proper position.

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